

Peroxiredoxin 4 (PRDX4) Rabbit Polyclonal Antibody

Information

Catalog Number: ZH4089

Host: Rabbit

Reactivity: Human;Mouse;Rat

Applications: WB;IHC-P;IF

Recommended Dilutions:

WB:1:500-1:2000

IHC-P:1:100-1:500

IF:1:50-1:200

UniProt ID:Q13162 O08807 Q9Z0V5

Predicted molecular weight: 31 kDa

Observed molecular weight:

Concentration:Please refer to the concentration information on the antibody vial label

Storage: Store at -20°C. Stable for one year after shipment.

Storage Notes: Aliquoting is unnecessary for -20°C storage.

Purity: Affinity-chromatography

Safety Notes: This product is for research use only, not for diagnostic or therapeutic use.